



CADD Prizm Pump Screen Message Glossary

Screen Message	Action
9-volt battery depleted/ install good battery	Replace battery. The pump will not start with a depleted battery.
9-volt battery low	The battery is low but the pump is operable. Change the battery soon.
Air in line detected	Air has been detected in the tubing; the tubing may contain air bubbles, or the tubing may not be threaded through the air detector. Press the “NEXT” key to silence the alarm, then: <ul style="list-style-type: none"> · Make sure the tubing is threaded properly. · If the tubing contains air bubbles, close the clamp on the IV tubing and disconnect the tubing from your IV. Remove the tubing from the pump and “flick” the tubing with your finger to break up the bubbles.
Cassette not attached / pump will not run	The pump will not start without a cassette attached. Attach the cassette and start the pump as directed.
Error detected/E (code)	A pump fault has occurred. Close the tubing clamp and turn off the pump. Call your pharmacy.
High pressure	The pump has detected high pressure which may be caused by a downstream blockage, kink in the tubing or a closed clamp on the tubing or IV. <ul style="list-style-type: none"> · Press the “Start/Stop” key to stop the pump and silence the alarm. Remove the obstruction and restart the pump by pressing the “Start/Stop” key.
Reservoir volume is zero	The reservoir volume has reached 0.0 mL which means that the medication has completely infused. Prepare to start a new medication bag as instructed.
Reservoir volume is low	The reservoir volume is low meaning it is almost time to change your medication bag. Prepare to install a new medication bag as instructed.